
Instructions for use

RAT CONTROL SET

REF**AR K-8000R**

96

RUO

For research
use only –
Not for use
in diagnostic
procedures

RAT CONTROL SET

English

1. MATERIALS SUPPLIED, STORAGE AND STABILITY

Components	Content	Preparation	Storage at	Shelf life
Rat Control 1/2: Two levels of Testosterone, Corticosterone, Prolactin, TSH, Progesterone and Estradiol in rat serum	1 Set 2 vials à 1 ml	lyophilized; reconstitute each vial with 1 ml dist. H ₂ O 30 minutes before use	2 - 8 °C; after reconstitution store in aliquots at -20 °C	6 months after reconstitution or until the expiration date

2. APPLICATION

According to the respective package insert of the ELISA kits

AR E-8000 Testosterone rat/mouse

AR E-8100 Corticosterone rat/mouse

AR E-8200 Prolactin rat

AR E-8600 TSH rat

AR E-8700 Progesterone rat/mouse

AR E-8800 Estradiol rat

the control can be used for internal quality control.

The procedure calls for **equivalent sample volume** per test.

3. TARGET VALUES

The values provided here are intended as guidelines. Because of many sources of inter-laboratory variations each laboratory should develop its own expected values.

Parameter		Mean (ng/ml)	Target Range (ng/ml)
Testosterone	Control 1	2.9	1.7 - 4.0
	Control 2	6.5	3.9 - 9.2

Parameter		Mean (ng/ml)	Target Range (ng/ml)
Corticosterone	Control 1	45.3	27.2 - 63.5
	Control 2	214	128 - 300

Parameter		Mean (ng/ml)	Target Range (ng/ml)
Prolactin	Control 1	15.8	9.5 - 22.1
	Control 2	27.2	16.3 - 38.1


Parameter		Mean (ng/ml)	Target Range (ng/ml)
TSH	Control 1	2.6	1.6 - 3.6
	Control 2	9.9	6.0 - 13.9

Parameter		Mean (ng/ml)	Target Range (ng/ml)
Progesterone	Control 1	6.0	3.5 - 8.4
	Control 2	17.9	10.7 - 25.0

Parameter		Mean (pg/ml)	Target Range (pg/ml)
Estradiol	Control 1	110	65 - 154
	Control 2	241	145 - 337

 **For updated literature or any other information please contact your local supplier.**

Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date	LOT	Batch code		
	Consult instructions for use	CONT	Content		
	Caution	REF	Catalogue number	RUO	For research use only!

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