


**Instructions for use**  
**Glutamate ELISA**

Please use only the valid version of the Instructions for Use provided with the kit

**REF****BA E-2400R**

96

**RUO**

For research  
use only –  
Not for use  
in diagnostic  
procedures

## Glutamate ELISA

### 1. Introduction

#### 1.1 Intended use and principle of the test

Enzyme Immunoassay for the quantitative determination of L-glutamate in urine and various biological samples.

After extraction and derivatisation glutamate is quantitatively determined by ELISA.

The competitive ELISA uses the microtiter plate format. The antigen is bound to the solid phase of the microtiter plate. The derivatized analyte concentrations in the standards, controls and samples and the solid phase bound analyte compete for a fixed number of antibody binding sites. When the system is in equilibrium, free antigen and free antigen-antibody complexes are removed by washing. The antibody bound to the solid phase is detected by an anti-rabbit IgG-peroxidase conjugate using TMB as a substrate. The reaction is monitored at 450 nm.

Quantification of unknown samples is achieved by comparing their absorbance with a standard curve prepared with known standards.

### 2. Procedural cautions, guidelines, warnings and limitations

#### 2.1 Procedural cautions, guidelines and warnings

- (1) This kit is intended for research use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- (2) The principles of Good Laboratory Practice (GLP) have to be followed.
- (3) In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
- (4) All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
- (5) For dilution or reconstitution purposes, use deionized, distilled, or ultra-pure water.
- (6) The microplate contains snap-off strips. Unused wells must be stored at 2 – 8 °C in the sealed foil pouch with desiccant and used in the frame provided.
- (7) Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
- (8) Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
- (9) Incubation times do influence the results. All wells should be handled in the same order and time intervals.
- (10) To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
- (11) A standard curve must be established for each run.
- (12) The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report provided with the kit.
- (13) Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.
- (14) Avoid contact with Stop Solution containing 0.25 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns. In case of contact with eyes or skin, rinse off immediately with water.
- (15) TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them.
- (16) For information on hazardous substances included in the kit please refer to Safety Data Sheet (SDS). The Safety Data Sheet for this product is made available directly on the website of the manufacturer or upon request.
- (17) Kit reagents must be regarded as hazardous waste and disposed of according to national regulations.
- (18) In case of any severe damage to the test kit or components, the manufacturer has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components must not be used for a test run. They must be stored properly until the manufacturer decides what to do with them. If it is decided that they are no longer suitable for measurements, they must be disposed of in accordance with national regulations.

#### 2.2 Limitations

Any inappropriate handling of samples or modification of this test might influence the results.

##### 2.2.1 Interfering substances

Avoid excess of acid: excess of acid might exceed the buffer capacity of the dilution buffer. A **pH of 5.0** during the extraction is mandatory.

### 2.2.2 Drug interferences

There are no known substances (drugs, food) which ingestion interferes with the measurement of glutamate level in the sample.

### 2.2.3 High-Dose-Hook effect


No hook effect was observed in this test.

## 3. Storage and stability

Store the unopened reagents at 2 – 8 °C until expiration date. Do not use components beyond the expiry date indicated on the kit labels. Once opened the reagents are stable for 2 months when stored at 2 – 8 °C. Once the resealable pouch has been opened, care should be taken to close it tightly with desiccant again.

## 4. Materials

### 4.1 Contents of the kit

<b>BA D-0090</b>	<b>FOILS</b>	<b>Adhesive Foil</b> – Ready to use
Contents:	Adhesive Foils in a resealable pouch	
Volume:	1 x 4 foils	
<b>BA D-0024</b>	<b>REAC-PLATE</b>	<b>Reaction Plate</b> – Ready to use
Contents:	1 x 96 well plate, empty in a resealable pouch	
<b>BA E-2442</b>	<b>EXTRACT-PLATE 48</b>	<b>Extraction Plate</b> – Ready to use
Contents:	2 x 48 well plate, precoated with cation exchanger in a resealable pouch	
<b>BA E-0030</b>	<b>WASH-CONC 50x</b>	<b>Wash Buffer Concentrate</b> – Concentrated 50x
Contents:	Buffer with a non-ionic detergent and physiological pH	
Volume:	1 x 20 ml/vial, light purple cap	
<b>BA E-0040</b>	<b>CONJUGATE</b>	<b>Enzyme Conjugate</b> – Ready to use
Contents:	Goat anti-rabbit immunoglobulins conjugated with peroxidase	
Volume:	1 x 12 ml/vial, red cap	
<b>BA E-0055</b>	<b>SUBSTRATE</b>	<b>Substrate</b> – Ready to use
Contents:	Chromogenic substrate containing tetramethylbenzidine, substrate buffer and hydrogen peroxide	
Volume:	1 x 12 ml/black vial, black cap	
<b>BA E-0080</b>	<b>STOP-SOLN</b>	<b>Stop Solution</b> – Ready to use
Contents:	0.25 M sulfuric acid	
Volume:	1 x 12 ml/vial, light grey cap	
Hazards identification:	 H290 May be corrosive to metals.	
<b>BA E-2431</b>	<b>GLUT</b>	<b>Glutamate Microtiter Strips</b> – Ready to use
Contents:	1 x 96 well (12x8) antigen precoated microwell plate in a resealable foil pouch with desiccant	
<b>BA E-2410</b>	<b>AS GLUT</b>	<b>Glutamate Antiserum</b> – Ready to use
Contents:	Rabbit anti-glutamate antibody, blue coloured	
Volume:	1 x 6 ml/vial, blue cap	
<b>BA E-2413</b>	<b>ASSAY-BUFF</b>	<b>Assay Buffer</b> – Ready to use
Contents:	Buffer with alkaline pH	
Volume:	1 x 20 ml/vial, yellow cap	

**BA E-2428** **EQUA-REAG** **Equalizing Reagent** – Lyophilized

Contents: Lyophilized protein

Volume: 1 vial, brown cap

**Standards and Controls** – Ready to use

Cat. no.	Component	Colour/Cap	Concentration µg/ml	Concentration µmol/l	Volume/ Vial
<b>BA E-2401</b>	<b>STANDARD A</b>	white	0	0	4 ml
<b>BA E-2402</b>	<b>STANDARD B</b>	light yellow	0.6	4.08	4 ml
<b>BA E-2403</b>	<b>STANDARD C</b>	orange	2	13.6	4 ml
<b>BA E-2404</b>	<b>STANDARD D</b>	dark blue	6	40.8	4 ml
<b>BA E-2405</b>	<b>STANDARD E</b>	light grey	20	136	4 ml
<b>BA E-2406</b>	<b>STANDARD F</b>	black	60	408	4 ml
<b>BA E-2451</b>	<b>CONTROL 1</b>	light green	Refer to QC-Report for expected value and acceptable range!		4 ml
<b>BA E-2452</b>	<b>CONTROL 2</b>	dark red			4 ml

Conversion: Glutamate (µg/ml) x 6.8 = Glutamate (µmol/l)

Contents: Acidic buffer with non-mercury preservative, spiked with defined quantity of Glutamate

**BA E-2446** **D-REAGENT** **D-Reagent** – Ready to use

Contents: Crosslinking agent in dimethylsulfoxide

Volume: 1 x 3 ml/vial, white cap

Hazards  
identification:

H317 May cause an allergic skin reaction.

**BA E-2458** **Q-BUFFER** **Q-Buffer** – Ready to use

Contents: TRIS buffer

Volume: 1 x 20 ml/vial, white cap

**BA E-2460** **DILUENT** **Diluent** – Ready to use

Contents: Buffer with sodium acetate

Volume: 1 x 20 ml/vial, dark green cap

**BA E-2787** **NAOH** **NaOH** – Ready to use

Contents: Sodium hydroxide solution

Volume: 1 x 2 ml/vial, purple cap

Hazards  
identification:

H290 May be corrosive to metals.

H315 Causes skin irritation.

H319 Causes serious eye irritation.

**4.2 Additional materials and equipment required but not provided in the kit**

- Calibrated precision pipettes to dispense volumes between 10 – 100 µl; 12.5 ml
- Microtiter plate shaker (shaking amplitude 3 mm; approx. 600 rpm)
- ELISA reader capable of reading absorbance at 450 nm and if possible 620 – 650 nm
- Shaker (shaking amplitude 3 mm; approx. 600 rpm)
- Absorbent material (paper towel)
- Vortex mixer
- Water (deionized, distilled, or ultra-pure)

**5. Sample collection and storage**

Various biological samples can be used for L-Glutamate determination. The assay was validated for urine samples.

## Urine

Urine stabilized with 10 µl 6 N HCl per 1 ml of urine sample can be used.

Storage: up to 6 hours (18 – 25 °C); up to 14 days (2 – 8 °C); up to 6 months (< -15 °C).

Repeated freezing and thawing should be avoided. Avoid exposure to direct sunlight.

## 6. Test procedure

Allow all reagents and samples to reach room temperature and mix thoroughly by gentle inversion before use. Duplicate determinations are recommended. It is recommended to number the strips of the microwell plate before usage to avoid any mix-up.

The binding of the antisera and of the enzyme conjugate and the activity of the enzyme are temperature dependent. The higher the temperature, the higher the absorption values will be. Varying incubation times will have similar influences on the absorbance. The optimal temperature during the Enzyme Immunoassay is between 20 – 25 °C.

During the overnight incubation at 2 – 8 °C with the antiserum, the temperature should be uniform all over the ELISA plate to avoid any drift and edge-effect.

⚠ In case of overflow, read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 405 nm.

### 6.1 Preparation of reagents

#### Wash Buffer

Dilute the 20 ml Wash Buffer Concentrate with water (deionized, distilled, or ultra-pure) to a final volume of 1000 ml.

Storage: 2 months at 2 – 8 °C

#### Equalizing Reagent

Reconstitute the Equalizing Reagent with **12.5 ml** of **Assay Buffer**.

Reconstituted Equalizing Reagent which is not used immediately has to be stored in aliquots for max 2 months at -20 °C and may be thawed only once.

#### D-Reagent

The D-Reagent has a freezing point of 18.5 °C. To ensure that the D-Reagent is liquid when being used, it must be ensured that the D-Reagent has reached room temperature and forms a homogeneous, crystal-free solution.

#### Glutamate Microtiter Strips

In rare cases residues of the blocking and stabilizing reagent can be seen in the wells as small, white dots or lines. These residues do not influence the quality of the product.

#### Extraction Plate

In rare cases residues of the cation exchanger can be seen in the wells as small, black dots or lines. These residues do not influence the quality of the product.

### 6.2 Preparation of samples

The Glutamate ELISA is a flexible test system for various biological sample types and volumes. It is not possible to give a general advice how to prepare the samples. However, the following basics should help the researcher to adapt the protocol to his specific needs:

- Avoid excess of acid: excess of acid might exceed the buffer capacity of the dilution buffer. A **pH of 5.0** during the extraction is mandatory.
- It is advisable to perform a **Proof of Principle** to determine the recovery of glutamate from the samples. Prepare a stock solution of glutamate. Add small amounts (to change the native sample matrix as less as possible) of the stock solutions to the sample matrix and check the recovery.
- The sample volume determines the sensitivity of this test. Determine the sample volume needed to determine glutamate in your sample by testing different amounts of sample volumes.
- If a sample volume **< 100 µl** is used, water (deionized, distilled, or ultra-pure) has to be added to a final **volume of 100 µl**.

*If you need any support in establishing a protocol for your specific purposes, do not hesitate to contact the manufacturer directly!*

## Extraction

1.	Pipette <b>100 µl</b> of the <b>standards, controls</b> and <b>samples</b> into the appropriate wells of the <b>Extraction Plate</b> .
2.	Add <b>100 µl</b> of the <b>Diluent</b> to all wells. Cover plate with <b>Adhesive Foil</b> and <b>shake</b> for <b>10 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
3.	Use <b>25 µl</b> for the subsequent <b>derivatization</b> !

## Derivatization

1.	Pipette <b>25 µl</b> of the <b>extracted standards, controls</b> and <b>samples</b> into the appropriate wells of the <b>Reaction Plate</b> .
2.	Pipette <b>10 µl</b> of <b>NaOH</b> into all wells.
3.	Pipette <b>50 µl</b> of the <b>Equalizing Reagent</b> into all wells.
4.	Pipette <b>10 µl</b> of the <b>D-Reagent</b> into all wells.
5.	Cover plate with <b>Adhesive Foil</b> and shake for <b>2 h</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
6.	Pipette <b>75 µl</b> of the <b>Q-Buffer</b> into all wells.
7.	Shake for <b>10 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
8.	Use <b>25 µl</b> for the <b>ELISA</b> !

## 6.3 Glutamate ELISA

1.	Pipette <b>25 µl</b> of the <b>prepared standards, controls and samples</b> into the appropriate wells of the <b>Glutamate Microtiter Strips</b> .
2.	Pipette <b>50 µl</b> of the <b>Glutamate Antiserum</b> into all wells and mix shortly.
3.	Cover plate with <b>Adhesive Foil</b> and incubate for <b>15 – 20 h</b> (overnight) at <b>2 – 8 °C</b> .
4.	Remove the foil. Discard or aspirate the content of the wells. Wash the plate <b>3 x</b> by adding <b>300 µl</b> of <b>Wash Buffer</b> , <b>discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.
5.	Pipette <b>100 µl</b> of the <b>Enzyme Conjugate</b> into all wells.
6.	Incubate for <b>30 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
7.	Discard or aspirate the contents of the wells and wash the plate <b>3 x</b> by adding <b>300 µl</b> of <b>Wash Buffer</b> , <b>discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.
8.	Pipette <b>100 µl</b> of the <b>Substrate</b> into all wells and incubate for <b>20 – 30 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm). <b>Avoid exposure to direct sunlight!</b>
9.	Add <b>100 µl</b> of the <b>Stop Solution</b> to each well and shake the microtiter plate to ensure a homogeneous distribution of the solution.
10.	<b>Read</b> the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to <b>450 nm</b> (if available a reference wavelength between 620 nm and 650 nm is recommended).

## 7. Calculation of results

Measuring range	Glutamate
	0.26 – 60 µg/ml

The standard curve is obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis). Use non-linear regression for curve fitting (e. g. 4-parameter, marquardt).

⚠ *This assay is a competitive assay. This means: the OD-values are decreasing with increasing concentrations of the analyte. OD-values found below the standard curve correspond to high concentrations of the analyte in the sample and have to be reported as being positive.*

## Samples and Controls

The concentrations of the samples (100 µl undiluted sample used) and controls can be read directly from the standard curve.

⚠ In case < 100 µl sample volume was used, concentrations of the samples taken from the standard curve have to be multiplied by a correction factor:

$$\text{Correction factor} = \frac{100 \mu\text{l (volume of standards)}}{\text{sample volume } (\mu\text{l})}$$

⚠ In case samples were pre-diluted correct the read values for the pre-dilution

## Conversion


Glutamate (µg/ml) x 6.8 = Glutamate (µmol/l)

## 7.1 Quality control

The confidence limits of the kit controls are indicated on the QC-Report.

## 7.2 Typical standard curve

⚠ Example, do not use for calculation!



## 8. Assay characteristics

Various biological samples can be used for L-Glutamate determination. The assay was validated for urine samples.

Analytical Sensitivity	Glutamate
Limit of Blank (LOB)	0.11 µg/ml
Limit of Detection (LOD)	0.17 µg/ml
Limit of Quantification (LOQ)	0.26 µg/ml

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)
		Glutamate
	L-Glutamine	< 0.4
	Glycine	< 0.4
	β-Alanine	< 0.4
	L-Alanine	< 0.4
	L-Aspartic Acid	< 0.4
	GABA	< 0.4
	5-Amino-n-valeric Acid	< 0.4













Precision							
Intra-Assay				Inter-Assay			
Sample	n	Mean $\pm$ SD ( $\mu\text{g/ml}$ )	CV (%)	Sample	n	Mean $\pm$ SD ( $\mu\text{g/ml}$ )	CV (%)
1	10	$0.8 \pm 0.1$	10.8	1	13	$1.7 \pm 0.24$	14.3
2	10	$1.3 \pm 0.1$	8.7	2	14	$5.0 \pm 0.57$	11.4
3	10	$2.2 \pm 0.1$	6.3	3	14	$10.6 \pm 0.73$	6.9
4	10	$4.8 \pm 0.2$	4.0	4	13	$3.0 \pm 0.43$	14.2
5	10	$12.5 \pm 0.6$	4.6	5	14	$5.6 \pm 0.71$	12.5
6	10	$39.7 \pm 2.2$	5.6	6	14	$10.0 \pm 0.87$	8.7

Linearity		Serial dilution up to	Range (%)	Mean (%)
	Urine	1:64	94 – 113	105

Recovery		Range ( $\mu\text{g/ml}$ )	Range (%)	Mean (%)
	Urine	1.25 – 41.0	97 – 108	102

 **For literature or any other information please contact your local supplier.**

#### Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Use-by date		Batch code		
	Consult instructions for use		Content		
	Caution		Catalogue number		Distributor
	Date of manufacture				For research use only!