
**Instructions for use**  
**GABA ELISA**

Please use only the valid version of the Instructions for Use provided with the kit

**REF****BA E-2500****IVD****CE**

## 1. Introduction

### 1.1 Intended use and principle of the test

Enzyme Immunoassay for the quantitative determination of gamma-aminobutyric acid (GABA) in urine, to evaluate GABA homeostasis. The determination of gamma-aminobutyric acid (GABA) in urine is helpful for the determination of neurostress.

After extraction and derivatization gamma-aminobutyric acid (GABA) is quantitatively determined by ELISA. The competitive ELISA uses the microtiter plate format. The antigen is bound to the solid phase of the microtiter plate. The derivatized analyte concentrations of the standards, controls and samples and the solid phase bound analyte compete for a fixed number of antiserum binding sites. When the system is in equilibrium, free antigen and free antigen-antibody complexes are removed by washing. The antibody bound to the solid phase is detected by an anti-rabbit IgG-peroxidase conjugate using TMB as a substrate. The reaction is monitored at 450 nm.

Quantification of unknown samples is achieved by comparing their absorbance with a standard curve prepared with known standards.

### 1.2 Clinical application

GABA ( $\gamma$ -aminobutyric acid) is one of the most important inhibitory neurotransmitter in the central nervous system (CNS). GABA operates through interneurons by the inhibition of the release of excitatory neurotransmitters. In the CNS GABA is synthesized from L-glutamic acid which is an excitatory neurotransmitter.

Many publications postulate, that the determination of GABA in urine is helpful for the determination of neurostress. The collective term „neurostress“ refers to many physical and psychological complaints caused due to our modern way of life, unfavourable environmental conditions, poor diet, medications, occupational and social stress, sleep deprivation, overstimulation or genetic predisposition. The resulting symptoms are burnout, depression, insomnia, chronic fatigue syndrome, fibromyalgia, multiple chemical sensitivity and other chronic pathologies. Furthermore, many publications describe an impaired or dysregulated glutamic acid balance in relation to previously listed symptoms. Based on this some laboratories offer the determination of specific neurostress profiles including, among others, the determination of GABA. In this case, GABA is detected primarily in the second morning urine by several methods. To generate a profound neurostress profile, also other analytes (glutamate, serotonin, norepinephrine, dopamine, epinephrine, melatonin, DHEA and cortisol) should be determined.

Determined GABA-levels beyond the reference value should be discussed and clarified with a therapist or an attending physician.

## 2. Procedural cautions, guidelines, warnings and limitations

### 2.1 Procedural cautions, guidelines and warnings

- (1) This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- (2) This assay was validated for a certain type of sample as indicated in *Intended Use* (please refer to Chapter 1). Any off-label use of this kit is in the responsibility of the user and the manufacturer cannot be held liable.
- (3) The principles of Good Laboratory Practice (GLP) have to be followed.
- (4) In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
- (5) All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
- (6) For dilution or reconstitution purposes, use deionized, distilled, or ultra-pure water.
- (7) The microplate contains snap-off strips. Unused wells must be stored at 2 – 8 °C in the sealed foil pouch with desiccant and used in the frame provided.
- (8) Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
- (9) Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
- (10) Incubation times do influence the results. All wells should be handled in the same order and time intervals.
- (11) To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
- (12) A standard curve must be established for each run.
- (13) The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report provided with the kit.
- (14) Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.

- (15) Avoid contact with Stop Solution containing 0.25 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns. In case of contact with eyes or skin, rinse off immediately with water.
- (16) TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them.
- (17) For information on hazardous substances included in the kit please refer to Safety Data Sheet (SDS). The Safety Data Sheet for this product is made available directly on the website of the manufacturer or upon request.
- (18) Kit reagents must be regarded as hazardous waste and disposed of according to national regulations.
- (19) The expected reference values reported in this test instruction are only indicative. It is recommended that each laboratory establishes its own reference intervals.
- (20) In case of any severe damage to the test kit or components, the manufacturer has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components must not be used for a test run. They must be stored properly until the manufacturer decides what to do with them. If it is decided that they are no longer suitable for measurements, they must be disposed of in accordance with national regulations.

**2.2 Limitations**

Any inappropriate handling of samples or modification of this test might influence the results.

**2.2.1 Interfering substances**

**Urine**

Please note the sample preparation stabilization of the urine sample! It cannot be excluded that high acid concentrations lead to incorrect results. Up to 20 µl 6 N HCl per 1 ml urine no influence on the results was observed.

**2.2.2 Drug interferences**

It is recommended to avoid food supplements which might influence GABA levels (lemon balm, valerian, vitamin B6, L-theanine and kava) 24 hours before urine sampling. Furthermore, additional diet restrictions before urine sampling are not described in the literature. In case of unusual GABA levels it is recommended to check if these are caused by interaction of the above listed substances.

**2.2.3 High-Dose-Hook effect**



No hook effect was observed in this test.

**3. Storage and stability**

Store the unopened reagents at 2 – 8 °C until expiration date. Do not use components beyond the expiry date indicated on the kit labels. Once opened the reagents are stable for 2 months when stored at 2 – 8 °C. Reconstituted Equalizing Reagent which is not used immediately has to be stored in aliquots for max. 2 months at < -15 °C and may be thawed only once. Once the resealable pouch has been opened, care should be taken to close it tightly with desiccant again.

**4. Materials**

**4.1 Contents of the kit**

- BA D-0090**  **Adhesive Foil** – Ready to use  
 Contents: Adhesive Foils in a resealable pouch  
 Volume: 3 x 4 foils
- BA D-0033**  **Macrotiter Plate** – Ready to use  
 Contents: 2 x 48 well plate, empty in a resealable pouch
- BA E-2442**  **Extraction Plate** – Ready to use  
 Contents: 2 x 48 well plate, precoated with cation exchanger in a resealable pouch
- BA E-0030**  **Wash Buffer Concentrate** – Concentrated 50x  
 Contents: Buffer with a non-ionic detergent and physiological pH  
 Volume: 1 x 20 ml/vial, light purple cap

**BA E-0040** **CONJUGATE** **Enzyme Conjugate** – Ready to use  
 Contents: Goat anti-rabbit immunoglobulins conjugated with peroxidase  
 Volume: 1 x 12 ml/vial, red cap

**BA E-0055** **SUBSTRATE** **Substrate** – Ready to use  
 Contents: Chromogenic substrate containing tetramethylbenzidine, substrate buffer and hydrogen peroxide  
 Volume: 1 x 12 ml/black vial, black cap

**BA E-0080** **STOP-SOLN** **Stop Solution** – Ready to use  
 Contents: 0.25 M sulfuric acid  
 Volume: 1 x 12 ml/vial, light grey cap

Hazards identification:   
 H290 May be corrosive to metals.

**BA E-2531** **GABA** **GABA Microtiter Strips** – Ready to use  
 Contents: 1 x 96 well (12x8) antigen precoated microwell plate in a resealable foil pouch with desiccant

**BA E-2510** **AS GABA** **GABA Antiserum** – Ready to use  
 Contents: Rabbit anti-GABA antibody, blue coloured  
 Volume: 1 x 6 ml/vial, blue cap

**Standards and Controls** – Ready to use

Cat. no.	Component	Colour/Cap	Concentration ng/ml	Concentration nmol/l	Volume/Vial
<b>BA E-2501</b>	<b>STANDARD A</b>	white	0	0	4 ml
<b>BA E-2502</b>	<b>STANDARD B</b>	light yellow	75	727	4 ml
<b>BA E-2503</b>	<b>STANDARD C</b>	orange	250	2,425	4 ml
<b>BA E-2504</b>	<b>STANDARD D</b>	dark blue	750	7,275	4 ml
<b>BA E-2505</b>	<b>STANDARD E</b>	light grey	2,500	24,250	4 ml
<b>BA E-2506</b>	<b>STANDARD F</b>	black	7,500	72,750	4 ml
<b>BA E-2551</b>	<b>CONTROL 1</b>	light green	Refer to QC-Report for expected value and acceptable range!		4 ml
<b>BA E-2552</b>	<b>CONTROL 2</b>	dark red			4 ml

Conversion: GABA (ng/ml) x 9.7 = GABA (nmol/l)

Contents: Acidic buffer with non-mercury preservative, spiked with defined quantity of GABA

**BA E-2413** **ASSAY-BUFF** **Assay Buffer** – Ready to use  
 Contents: Buffer with alkaline pH  
 Volume: 1 x 20 ml/vial, yellow cap

**BA E-2428** **EQUA-REAG** **Equalizing Reagent** – Lyophilized  
 Contents: Lyophilized protein  
 Volume: 1 vial, brown cap

**BA E-2446** **D-REAGENT** **D-Reagent** – Ready to use  
 Contents: Crosslinking agent in dimethylsulfoxide  
 Volume: 1 x 3 ml/vial, white cap

Hazards identification:   
 H317 May cause an allergic skin reaction.

**BA E-2458**    **Q-BUFFER**                    **Q-Buffer** – Ready to use

Contents:    TRIS buffer

Volume:     1 x 20 ml/vial, white cap

**BA E-2561**    **I-BUFFER**                    **I-Buffer** – Concentrated

Contents:    Buffer with non-ionic detergent and non-mercury preservative

Volume:     1 x 4 ml/vial, light red cap

**BA E-2541**    **ELUTION-BUFF**                **Elution-Buffer** – Ready to use

Contents:    Buffer with citric acid

Volume:     1 x 20 ml/vial, dark green cap

**BA E-2560**    **DILUENT**                    **Diluent** – Ready to use

Contents:    Buffer with acidic pH

Volume:     2 x 20 ml/vial, blue cap

**BA E-2787**    **NAOH**                        **NaOH** – Ready to use

Contents:    Sodium hydroxide solution

Volume:     1 x 2 ml/vial, purple cap

Hazards  
identification:

H290 May be corrosive to metals.  
H315 Causes skin irritation.  
H319 Causes serious eye irritation.

#### 4.2 Additional materials and equipment required but not provided in the kit

- Calibrated precision pipettes to dispense volumes between 10 – 500 µl; 12.5 ml
- Microtiter plate washing device (manual, semi-automated or automated)
- ELISA reader capable of reading absorbance at 450 nm and if possible 620 – 650 nm
- Microtiter plate shaker (shaking amplitude 3 mm, approx. 600 rpm)
- Vortex mixer
- Absorbent material (paper towel)
- Water (deionized, distilled, or ultra-pure)

#### 5. Sample collection and storage

##### Urine

Spontaneous urine (second morning urine) stabilized with 10 µl 6 N HCl per 1 ml of urine sample can be used. The measurement results are related to the creatinine content of the sample.

Storage: up to 6 hours (18 – 25 °C); up to 14 days (2 – 8 °C); up to 6 months (< -15 °C).

Repeated freezing and thawing should be avoided. Avoid exposure to direct sunlight.

#### 6. Test procedure

Allow all reagents and samples to reach room temperature and mix thoroughly by gentle inversion before use. Duplicate determinations are recommended. It is recommended to number the strips of the microwell plate before usage to avoid any mix-up.

The binding of the antisera and of the enzyme conjugate and the activity of the enzyme are temperature dependent. The higher the temperature, the higher the absorption values will be. Varying incubation times will have similar influences on the absorbance. The optimal temperature during the Enzyme Immunoassay is between 20 – 25 °C.

During the overnight incubation at 2 – 8 °C with the antiserum, the temperature should be uniform all over the ELISA plate to avoid any drift and edge-effect.

⚠ *The use of a microtiter plate shaker with the following specifications is mandatory: shaking amplitude 3 mm; approx. 600 rpm. Shaking with differing settings might influence the results.*

## 6.1 Preparation of reagents

### Wash Buffer

Dilute the 20 ml Wash Buffer Concentrate with water (deionized, distilled, or ultra-pure) to a final volume of 1000 ml.

Storage: 2 months at 2 – 8 °C

### Equalizing Reagent

Reconstitute the Equalizing Reagent with **12.5 ml** of **Assay Buffer**.

Reconstituted Equalizing Reagent which is not used immediately has to be stored in aliquots for max. 2 months at < -15 °C and may be thawed only once.

### I-Buffer

Dilute the 4 ml I-Buffer Concentrate with water (deionized, distilled, or ultra-pure) to a final volume of 400 ml.

Storage: 2 months at 2 – 8 °C

### D-Reagent

The D-Reagent has a freezing point of 18.5 °C. Make sure that the D-Reagent has reached room temperature and forms a homogeneous, crystal-free solution.

### GABA Microtiter Strips

In rare cases residues of the blocking and stabilizing reagent can be seen in the wells as small, white dots or lines. These residues do not influence the quality of the product.

### Extraction Plate

In rare cases residues of the cation exchanger can be seen in the wells as small, black dots or lines. These residues do not influence the quality of the product.

## 6.2 Preparation of samples

### Extraction

<b>1.</b> Pipette <b>100 µl</b> of the <b>standards, controls</b> and <b>samples</b> into the appropriate wells of the <b>Extraction Plate</b> .
<b>2.</b> Add <b>100 µl</b> of the <b>Diluent</b> to all wells. Cover plate with <b>Adhesive Foil</b> and incubate for <b>15 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
<b>3.</b> <b>Discard</b> and blot dry by tapping the inverted plate on absorbent material. <b>Add 500 µl</b> of <b>I-Buffer</b> to each well and incubate for <b>5 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
<b>4.</b> <b>Discard</b> the wash and blot dry by tapping the inverted plate on absorbent material.
<b>5.</b> Pipette <b>150 µl</b> of <b>Elution Buffer</b> into the appropriate wells of the <b>Extraction Plate</b> . Cover plate with <b>Adhesive Foil</b> and incubate for <b>10 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
<b>6.</b> Use <b>100 µl</b> for the subsequent <b>derivatization!</b>

### Derivatization

<b>1.</b> Pipette <b>100 µl</b> of the <b>extracted standards, controls</b> and <b>samples</b> into the appropriate wells of the <b>Macrotiter Plate</b> .
<b>2.</b> Pipette <b>10 µl</b> of the <b>NaOH</b> into all wells.
<b>3.</b> Add <b>50 µl</b> of the <b>Equalizing Reagent</b> (fresh prepared before assay) to all wells.
<b>4.</b> Incubate for <b>1 min</b> on a shaker (approx. 600 rpm).
<b>5.</b> Pipette <b>10 µl</b> of the <b>D-Reagent</b> into all wells.
<b>6.</b> Cover plate with <b>Adhesive Foil</b> and incubate for <b>2 h</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
<b>7.</b> Pipette <b>100 µl</b> <b>Q-Buffer</b> into all wells.
<b>8.</b> Shake for <b>10 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
<b>9.</b> Use <b>50 µl</b> for the subsequent <b>ELISA!</b>

**GABA ELISA**

<b>1.</b>	Pipette <b>50 µl</b> of the <b>derivatized standards, controls</b> and <b>samples</b> into the appropriate wells of the <b>GABA Microtiter Strips</b> .
<b>2.</b>	Pipette <b>50 µl</b> of the <b>GABA Antiserum</b> into all wells and mix shortly.
<b>3.</b>	Cover plate with <b>Adhesive Foil</b> and incubate for <b>15 – 20 h</b> (overnight) at <b>2 – 8 °C</b> .
<b>4.</b>	Remove the foil. Discard or aspirate the content of the wells. Wash the plate <b>3 x</b> by adding <b>300 µl</b> of <b>Wash Buffer, discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.
<b>5.</b>	Pipette <b>100 µl</b> of the <b>Enzyme Conjugate</b> into all wells.
<b>6.</b>	Incubate for <b>30 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
<b>7.</b>	Discard or aspirate the content of the wells. Wash the plate <b>3 x</b> by adding <b>300 µl</b> of <b>Wash Buffer, discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.
<b>8.</b>	Pipette <b>100 µl</b> of the <b>Substrate</b> into all wells and incubate for <b>20 – 30 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm). <b>Avoid exposure to direct sunlight!</b>
<b>9.</b>	Add <b>100 µl</b> of the <b>Stop Solution</b> to each well and shake the microtiter plate to ensure a homogeneous distribution of the solution.
<b>10.</b>	<b>Read</b> the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to <b>450 nm</b> (if available a reference wavelength between 620 nm and 650 nm is recommended).

**7. Calculation of results**

<b>Measuring range</b>	<b>GABA</b>
	59.1 – 7,500 ng/ml

The standard curve is obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis). Use non-linear regression for curve fitting (e.g. 4-parameter, marquardt).

⚠ *This assay is a competitive assay. This means: the OD-values are decreasing with increasing concentrations of the analyte. OD-values found below the standard curve correspond to high concentrations of the analyte in the sample and have to be reported as being positive.*

**Urine samples and controls**

The concentrations of the samples and controls can be read directly from the standard curve. Samples found with concentrations higher than the highest standard (Standard F) should be diluted accordingly with water (deionized, distilled, or ultra-pure) and have to be re-assayed.

**Conversion**

$GABA (ng/ml) \times 9.7 = GABA (nmol/l)$

**Expected reference value**

It is strongly recommended that each laboratory should determine its own reference value.

<b>Expected reference value</b>	<b>Urine</b>
	206 – 1,548 µg/g creatinine
	2 – 15 µmol/g creatinine
	0.23 – 1.7 mmol/mol creatinine


Values significantly outside the reference range should be assessed by a doctor.

**7.1 Quality control**

The confidence limits of the kit controls are indicated on the QC-Report.

## 7.2 Typical standard curve

⚠ Example, do not use for calculation!

## 8. Assay characteristics

Analytical Sensitivity	GABA
Limit of Blank (LOB)	19.6 ng/ml
Limit of Detection (LOD)	30.4 ng/ml
Limit of Quantification (LOQ)	59.1 ng/ml

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)
		GABA
	3-Aminobutanoic acid	< 0.1
	L-(+)-2-Aminobutyric acid	< 0.1
	β-Alanine	0.8
	L-Aspartic acid	< 0.1
	(S)-(+)-Glutamine	< 0.1
	Glycine	< 0.1
	L-Glutamic acid	< 0.1

Precision					
Intra-Assay			Inter-Assay		
Sample	Mean ± SD (ng/ml)	CV (%)	Sample	Mean ± SD (ng/ml)	CV (%)
1	135 ± 32	23.7	1	205 ± 25	12.4
2	241 ± 30	12.6	2	426 ± 40	9.4
3	451 ± 41	9.2	3	944 ± 69	7.3
4	1,021 ± 56	5.5	4	2,763 ± 195	7.1
5	2,871 ± 200	7.2			
6	6,327 ± 322	5.1			

Recovery		Range (ng/ml)	Range (%)	Mean (%)
	Urine		135 – 6,327	103 – 117

Linearity		Serial dilution up to	Range (%)	Mean (%)
	Urine		1:64	98 – 123

## 9. References/Literature

- (1) Bieger, W.P., *NeuroStress Guide*. 2011.
- (2) Bustillo, J.R., *Use of proton magnetic resonance spectroscopy in the treatment of psychiatric disorders: a critical update*. *Dialogues Clin Neurosci*, 2013. **15**(3): p. 329-37.
- (3) Duman, R.S., G. Sanacora, and J.H. Krystal, *Altered Connectivity in Depression: GABA and Glutamate Neurotransmitter Deficits and Reversal by Novel Treatments*. *Neuron*, 2019. **102**(1): p. 75-90.
- (4) Femenia, T., et al., *Dysfunctional hippocampal activity affects emotion and cognition in mood disorders*. *Brain Res*, 2012. **1476**: p. 58-70.
- (5) Flasnoecker, M., *Reise aus dem Stress - Körper, Geist und Psyche stärken*. 2015: W. Zuckschwerdt Verlag GmbH.
- (6) Frisardi, V., F. Panza, and A.A. Farooqui, *Late-life depression and Alzheimer's disease: the glutamatergic system inside of this mirror relationship*. *Brain Res Rev*, 2011. **67**(1-2): p. 344-55.
- (7) Gao, S.F. and A.M. Bao, *Corticotropin-releasing hormone, glutamate, and gamma-aminobutyric acid in depression*. *Neuroscientist*, 2011. **17**(1): p. 124-44.
- (8) Harris, R.E. and D.J. Clauw, *Imaging central neurochemical alterations in chronic pain with proton magnetic resonance spectroscopy*. *Neurosci Lett*, 2012. **520**(2): p. 192-6.
- (9) Kendell, S.F., J.H. Krystal, and G. Sanacora, *GABA and glutamate systems as therapeutic targets in depression and mood disorders*. *Expert Opin Ther Targets*, 2005. **9**(1): p. 153-68.
- (10) Krystal, J.H., et al., *Glutamate and GABA systems as targets for novel antidepressant and mood-stabilizing treatments*. *Mol Psychiatry*, 2002. **7 Suppl 1**: p. S71-80.
- (11) Lener, M.S., et al., *Glutamate and Gamma-Aminobutyric Acid Systems in the Pathophysiology of Major Depression and Antidepressant Response to Ketamine*. *Biol Psychiatry*, 2017. **81**(10): p. 886-897.
- (12) Sanacora, G., G. Treccani, and M. Popoli, *Towards a glutamate hypothesis of depression: an emerging frontier of neuropsychopharmacology for mood disorders*. *Neuropharmacology*, 2012. **62**(1): p. 63-77.
- (13) Strienz, J., *Nebennierenunterfunktion - Stress stört die Hormonbalance*. 3 ed. 2019: W. Zuckschwerdtverlag München.

For updated literature or any other information please contact your local supplier.

### Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Use-by date		Batch code		For in-vitro diagnostic use only!
	Consult instructions for use		Content		CE marking of conformity
	Caution		Catalogue number		Distributor
	Date of manufacture				

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
## 9. Referenzen/Literatur

- (1) Bieger, W.P., *NeuroStress Guide*. 2011.
- (2) Bustillo, J.R., *Use of proton magnetic resonance spectroscopy in the treatment of psychiatric disorders: a critical update*. Dialogues Clin Neurosci, 2013. **15**(3): p. 329-37.
- (3) Duman, R.S., G. Sanacora, and J.H. Krystal, *Altered Connectivity in Depression: GABA and Glutamate Neurotransmitter Deficits and Reversal by Novel Treatments*. Neuron, 2019. **102**(1): p. 75-90.
- (4) Femenia, T., et al., *Dysfunctional hippocampal activity affects emotion and cognition in mood disorders*. Brain Res, 2012. **1476**: p. 58-70.
- (5) Flasnoecker, M., *Reise aus dem Stress - Körper, Geist und Psyche stärken*. 2015: W. Zuckschwerdt Verlag GmbH.
- (6) Frisardi, V., F. Panza, and A.A. Farooqui, *Late-life depression and Alzheimer's disease: the glutamatergic system inside of this mirror relationship*. Brain Res Rev, 2011. **67**(1-2): p. 344-55.
- (7) Gao, S.F. and A.M. Bao, *Corticotropin-releasing hormone, glutamate, and gamma-aminobutyric acid in depression*. Neuroscientist, 2011. **17**(1): p. 124-44.
- (8) Harris, R.E. and D.J. Clauw, *Imaging central neurochemical alterations in chronic pain with proton magnetic resonance spectroscopy*. Neurosci Lett, 2012. **520**(2): p. 192-6.
- (9) Kendell, S.F., J.H. Krystal, and G. Sanacora, *GABA and glutamate systems as therapeutic targets in depression and mood disorders*. Expert Opin Ther Targets, 2005. **9**(1): p. 153-68.
- (10) Krystal, J.H., et al., *Glutamate and GABA systems as targets for novel antidepressant and mood-stabilizing treatments*. Mol Psychiatry, 2002. **7 Suppl 1**: p. S71-80.
- (11) Lener, M.S., et al., *Glutamate and Gamma-Aminobutyric Acid Systems in the Pathophysiology of Major Depression and Antidepressant Response to Ketamine*. Biol Psychiatry, 2017. **81**(10): p. 886-897.
- (12) Sanacora, G., G. Treccani, and M. Popoli, *Towards a glutamate hypothesis of depression: an emerging frontier of neuropsychopharmacology for mood disorders*. Neuropharmacology, 2012. **62**(1): p. 63-77.
- (13) Strienz, J., *Nebennierenunterfunktion - Stress stört die Hormonbalance*. 3 ed. 2019: W. Zuckschwerdtverlag München.

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