
**Instructions for use**  
**SALI SET 500**

Please use only the valid version of the Instructions for Use provided with the kit

**REF****SA D-6500****IVD****CE**

## 1. Introduction

### 1.1 Intended use

Sampling tubes for collection of saliva samples to be used in the LDN Saliva Immunoassays.

## 2. Procedural cautions, guidelines and warnings

### 2.1 Precautions, guidelines and warnings

- (1) Only the test instruction provided with the respective kit is valid.
- (2) Do not use components beyond expiry date as shown on the labels.
- (3) Be sure to follow the sample preparation instructions for the analysis product you are using.
- (4) Do not use other containers without validation. In general, plastic materials tend to adsorb steroids from saliva. The material used for the SALI SETs does not adsorb steroids in significant amounts.
- (5) Do not reuse.
- (6) Shipment: Human saliva samples are classified as potentially infective and fall into Category B. They must be classified as UN 3373 "Biological substance, Category B."
- (7) The unused product can be disposed of in household waste or commercial waste like household waste. Used devices must be regarded as hazardous waste and disposed of in accordance with national regulations.

## 3. Storage and stability

Store the collection system at room temperature until expiration date. Do not use components beyond the expiry date indicated on the label.

## 4. Materials

### 4.1 Contents of the kit

**SA D-6557**    **STRAWS**                      **Straws** - Ready to use

Content:        1 x 500\*, Polypropylene

**SA D-6558**    **TUBES**                                      **Tubes** - Ready to use

Content:        1 x 500\*, ultrapure polypropylene tubes

**SA D-6559**    **STICKERS**                                      **Stickers** - Ready to use

Content:        1 x 500\*, self-adhesive

\* approximately

## 5. Sample collection

### 5.1 Sampling Notes

⚠ *The requirements for sample collection and the procedure for returning the collected samples to the doctor/laboratory must be passed on to the patient in an appropriate form if the sample collection is not carried out under professional supervision.*

- (1) The concentration of many analytes in saliva is highly dynamic. Therefore, saliva samples must be collected at the times specified by the therapist.
- (2) If you want to check the average concentration of steroid hormones in saliva you have to take at least 5 samples within 2 hours. Aliquots of these 5 samples are mixed prior to analysis. Such a sampling strategy will result in reproducible results.
- (3) Collect Saliva **fasting or earliest 30 minutes after eating or drinking.**

⚠ *Depending on the analyte, additional restrictions and conditions for sample collection may exist. Please always check the instructions for sample preparation of the following analytical procedure and pass them on to your patients.*

- (4) Rinsing of the oral cavity with water is recommended just prior to saliva collection (especially after eating).
- (5) Do not chew gum, brush teeth, use mouthwashes, or smoke 30 minutes prior to saliva specimen collection.
- (6) Do **not use even the slightest red saliva samples** (because of blood contamination).

## 5.2 Sampling Procedure

- (1) Label the TUBE with label sticker.
- (2) Release the saliva with help of the piece of straw into the TUBE.
- (3) It is sufficient to fill half of the TUBE with saliva.



## 6. Pre-test setup instructions

In order to get a clear liquid please deep-freeze the saliva samples prior to use. After complete freezing thaw again, mix and centrifuge. If necessary, repeat this procedure several times, until the supernatant is a clear liquid and easy to pipette.

Please use only the valid version of the Instructions for Use provided with the kit

 **For updated literature or any other information please contact your local supplier.**

### Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date	<b>LOT</b>	Batch code	<b>IVD</b>	For in-vitro diagnostic use only!
	Consult instructions for use	<b>CONT</b>	Content	<b>CE</b>	CE labelled
	Caution	<b>REF</b>	Catalogue number		Single use only!

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